

WARNING: MAY FORM EXPLOSIVE PEROXIDES

Store, handle and dispose of per LBNL Controls for Peroxide Formers.

Keep in tightly closed original container. Avoid exposure to light, air and heat. If crystals, discoloration, or layering are visible, do not open.

Contact your EHS Health and Safety Rep for guidance.

Date Received_____ Date Opened_____ Testing Interval (months)_____

PEROXIDE TEST RESULTS

(If within the concentration range from 30-100 ppm material can be used but don't evaporate or concentrate)

Date_____ Result_____ Date_____ Result_____

Date_____ Result_____ Date_____ Result_____

Date_____ Result_____ Date_____ Result_____